


Coincidence summing correction in Genie-2000


Henrik Jäderström
Canberra Solutions AB
GammaSem-2010

- ▶ Introduction
- ▶ Coincidence summing algorithm
- ▶ Geometry description with LabSOCS
- ▶ Validation
- ▶ Summary



Introduction

- ▶ Genie-2000 is Canberra's software for gamma spectroscopy
- ▶ Coincidence summing correction first introduced in 2001
- ▶ Major upgrade in 2009, included summation with X-rays and annihilation photons
- ▶ Database updated and expanded



Algorithm

▶ **Coincidence summing correction factor (COI)**

$$COI_A = (1 - L^{\gamma-\gamma}_A - L^{\gamma-X,511}_A)(1 + S^{\gamma-\gamma}_A)$$

- ▶ $L^{\gamma-\gamma}_A$ and $S^{\gamma-\gamma}_A$ is the summing out and summing in probability for γ - γ coincidences
- ▶ $L^{\gamma-X,511}_A$ is the summing out probability for γ and EC X-ray and annihilation photons
- ▶ L and S are the sum of the partial decay chains for the gamma ray of interest

Summing in probability

- ▶ Summing in probability for a gamma chain with 3 gammas $A = Y_1 + Y_2 + Y_3$

$$S_A = \frac{\gamma_1}{\gamma_A} a_2 c_2 a_3 c_3 \frac{\epsilon_{p,1} \epsilon_{p,2} \epsilon_{p,3}}{\epsilon_{p,A}}$$

- ▶ γ is the absolute gamma yield, the number of photons of emitted per decay for the energy of interest.
- ▶ a is the branching ratio, the probability of particular transition from given nuclear state.
- ▶ c is the internal conversion factor, defined as the $1/(1+\alpha)$ where α is the ratio between the number of gammas and electrons emitted from a nuclear state.
- ▶ Summing in from more than 3 gammas can be neglected

Summing out probability

- ▶ More complicated than summing in
- ▶ Depends on the position of the transition of the decay chain
- ▶ Genie-2000 uses maximum 5 transitions in cascade
- ▶ For a cascade with 5 transitions and the transition of interest is the third the correction is, as given by de Corte

$$\begin{aligned} & \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot c_1 \cdot c_2 \cdot \Sigma_0 + \gamma_1/\gamma_2 \cdot a_2 \cdot c_2 \cdot \Sigma_1 + a_3 \cdot c_3 \cdot \Sigma_3 + a_3 \cdot a_4 \cdot c_4 \cdot \Sigma_4 - \\ & \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot c_1 \cdot c_2 \cdot \Sigma_0 \cdot \Sigma_2 - \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot a_3 \cdot c_2 \cdot c_3 \cdot \Sigma_0 \cdot \Sigma_3 - \gamma_0/\gamma_1 \cdot a_1 \cdot a_2 \cdot a_3 \cdot a_4 \cdot c_2 \cdot c_4 \cdot \Sigma_0 \cdot \Sigma_4 - \\ & \gamma_1/\gamma_2 \cdot a_2 \cdot a_3 \cdot c_2 \cdot c_3 \cdot \Sigma_1 \cdot \Sigma_3 - \gamma_1/\gamma_2 \cdot a_2 \cdot a_3 \cdot a_4 \cdot c_2 \cdot c_4 \cdot \Sigma_1 \cdot \Sigma_4 - a_3 \cdot a_4 \cdot c_3 \cdot c_4 \cdot \Sigma_3 \cdot \Sigma_4 + \\ & \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot a_3 \cdot c_1 \cdot c_2 \cdot c_3 \cdot \Sigma_0 \cdot \Sigma_1 \cdot \Sigma_3 + \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot a_3 \cdot a_4 \cdot c_1 \cdot c_2 \cdot c_4 \cdot \Sigma_0 \cdot \Sigma_1 \cdot \Sigma_4 + \\ & \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot a_3 \cdot a_4 \cdot c_2 \cdot c_3 \cdot c_4 \cdot \Sigma_0 \cdot \Sigma_3 \cdot \Sigma_4 + \gamma_1/\gamma_2 \cdot a_2 \cdot a_3 \cdot a_4 \cdot c_2 \cdot c_3 \cdot c_4 \cdot \Sigma_1 \cdot \Sigma_3 \cdot \Sigma_4 - \\ & \gamma_0/\gamma_2 \cdot a_1 \cdot a_2 \cdot a_3 \cdot a_4 \cdot c_1 \cdot c_2 \cdot c_3 \cdot c_4 \cdot \Sigma_0 \cdot \Sigma_1 \cdot \Sigma_3 \cdot \Sigma_4 \end{aligned} \quad (c)$$

- ▶ The total efficiency can be expressed as to account for coincidences with IC x-rays

IC and EC x-rays

- ▶ The total efficiency can be expressed as to account for coincidences with IC x-rays

$$\sum_i = \varepsilon_{t,i} + \alpha_K \omega_K \sum_j k_j \varepsilon_{t,k_j}$$

- ▶ α_K is the K-electron internal conversion factor, ω_K is the K-shell fluorescence yield, k_j is the KX-ray relative intensity
- ▶ For EC X-rays substitute in de Cortes formula

$$a_1 a_2 c_2 \frac{\gamma_0}{\gamma_2} \rightarrow \frac{\beta_0 a_0}{\beta_0 a_0 + \beta_1} = \frac{\beta_0 a_0}{T} = B_0, \quad \text{and} \quad \Sigma_0 \rightarrow P_{K,0} \omega_K \sum_i k_i \varepsilon_{t,k_i}$$


- ▶ $c = 1(1+\alpha_T)$ where α_T is the total internal conversion coefficient, β_i is the i-th decay branching ratio and $P_{K,j}$ is the K-electron capture probability of the EC decay

Annihilation photons

- ▶ Similar to EC x-ray correction but substitute

$$\beta_i \rightarrow \beta_i^{B+}, \quad \text{and} \quad P_{K,i} \omega_K \sum_j k_j \varepsilon_{t,k_j} \rightarrow 2 \times \varepsilon_{t,511},$$

- ▶ β_i^{B+} is the relative positron emission branching ratio and $\varepsilon_{i,511}$ is the total efficiency of the 511 keV photon and the 2 is for the two photons emitted from the positron annihilation



Subcascades

- ▶ Several nuclides contains more than one cascade and a partucular transition may be part of more than one of these subcascades
- ▶ 244 keV + 122 keV decay chain from Eu-152 is a part of seven subcascades
- ▶ In order not to overestimate the summing out probability it is necessary to reduce the total summing out probability for six subcascades involving this transition
- ▶ The total summing out probability can be expressed as

$$L_{244} = \sum_{i=1}^7 L_i - 6L_{244 \rightarrow 122}$$

Voluminous sources

- ▶ Activity calculated as $A=N_p/(\epsilon\gamma)$
- ▶ Dividing the source into n subsources gives the count rate from each subsource as $N_{p,i}=(N_p\epsilon_i)/(\epsilon n)$
- ▶ For a gamma line the coincidence corrected count rate is


$$N'_{p,g} = \frac{N_p}{\epsilon \cdot n} \sum_{i=1}^n \frac{\epsilon_i}{COI_i}$$

- ▶ The coincidence correction is then calculated as

$$COI = \frac{N_p}{N'_{p,g}} = \frac{\epsilon \cdot n}{\sum_{i=1}^n \frac{\epsilon_i}{COI_i}}$$

Total Efficiency


- ▶ In the calculations presented it is necessary to know the total efficiency for the gamma lines
- ▶ Genie offers two way of calculating the total efficiency
- ▶ Peak-To-Total calibrations (old way)
- ▶ Mathematical calculations with LabSOCS (preferred way)



Peak-To-Total Calibration


- ▶ **Peak Efficiency is normally already done**
- ▶ **Measure the number of counts in the full energy peak and the number of counts in the Compton continuum**
- ▶ **Usually extrapolated for the lowest energy**
- ▶ **Use mono-energetic sources**
- ▶ **From the interpolated curve and the peak efficiency the total efficiency can be calculated**

P/T calibration curve




Geometry description

- ▶ 3D geometry composer
- ▶ Efficiency calculations with LabSOCS/ISOCS algorithm
- ▶ 12 LabSOCS templates
- ▶ Customizable complex beaker template
- ▶ Energy range from 10 keV to 7000 keV
- ▶ Accuracy from 20% at 10 keV to 4 % at 2000 keV




Efficiency calculations

- ▶ Divide the geometry into 2^n voxels
- ▶ Place a point source in a random position in each voxel
- ▶ Calculate the efficiency for that point in space, including material attenuation
- ▶ Calculate efficiency for whole volume and all energies of interest
- ▶ Divide into 2^{n+1} voxels and redo calculations
- ▶ Repeat until convergence is reached
- ▶ Total efficiency is corrected for photon-buildup (scattering in attenuators) and albedo effects (scattering in surrounding materials)
- ▶ The build-up is based on MCNP simulations and is parameterized as a function of energy, material and geometry



Database used

- ▶ Data for the coincidence summing correction taken from ENSDF
- ▶ Atomic X-ray data
- ▶ Nuclides and gamma lines that are free from or subject to coincidence summing and their decay parameters
- ▶ Energies of gamma lines in summing-out chains
- ▶ Energies of gamma lines in summing-in chains
- ▶ Energies of gamma lines in sub-cascade chains
- ▶ Summarized in one database file containing 205 nuclides




J.K. Tuli, BNL-NCS051655-01/02-Rev

Validation


- ▶ **Extensive validation of the methods with different sources (all NIST traceable 1-2% uncertainty) and geometries**
- ▶ **Co-60, Y-88 and Ce-139 in three geometries: Close filter paper, 20 ml liquid vial and 2.8 l Marinelli beaker**
- ▶ **Ba-133, Eu-152 and Na-22 all disc sources**
- ▶ **All 205 nuclides and more than 2000 lines contained in the coincidence library compared to MCNP calculations**

- ▶ **Top: disc Na-22, Ba-133 and Eu-152 source**
- ▶ **Bottom: common mixed gamma calibration source with Ce-139, Y-88 and Co-60 shown**
- ▶ **Result**
 - ◆ **Gamma OK for Co-60 since no X-rays**
 - ◆ **Y-88, Ba-133 and Eu-152 are all better with X+g than with gamma only**



X-ray Cascade Summing Correction of Gammas v2.2

- ▶ Data are for Close Filter Paper (59 nuclides)
- ▶ Analysis using ISOCS TE for BE5030 [wide energy range detector]
- ▶ MCNP-CP used as reference
- ▶ Top – no CSC
- ▶ Middle – gamma only CSC
- ▶ Bottom – gamma & X-ray CSC
- ▶ Results: no bias, 2%sd slightly larger than for Gammas where there are no X-rays
- ▶ Similar results for 20 ml vial (15 nuclides) and Marinelli beaker (12 nuclides).
 - Mean: 0.5%
 - SD: 2.0%
 - %>±5%: 30%
 - %>±10%: 14%



Summary

- ▶ **Canberra has developed an algorithm to correct for coincidence summing in the software Genie-2000**
- ▶ **Can correct for gamma-gamma coincidences and gamma-x-ray/511 coincidences**
- ▶ **Does not require any sources for calibration**
- ▶ **Uses the LabSOCS code to calculate the both the peak and total efficiency (preferred way)**
- ▶ **Has undergone extensive validation which have shown that the systematic uncertainty from the correction is 5%**