

Mathematical efficiency calibration in gamma spectrometry for analysis of natural and fallout radionuclides in sediments

Daniela Pittauerová Bernd Hettwig
Helmut W. Fischer

University of Bremen, Institute for Environmental Physics, Germany

Sediment chronology



Efficiency calibration

Efficiency determination involves measurement of test sources and calculating efficiencies for given geometry and energy as

$$\epsilon_E = \frac{N_{net}}{t_L \cdot A \cdot I_\gamma},$$

where N_{net} is FEP net count, t_L is counting live time, A is source activity and I_γ probability of emission of the particular gamma line.

Efficiency calibration

Efficiency determination involves measurement of test sources and calculating efficiencies for given geometry and energy as

$$\epsilon_E = \frac{N_{net}}{t_L \cdot A \cdot I_\gamma},$$

where N_{net} is FEP net count, t_L is counting live time, A is source activity and I_γ probability of emission of the particular gamma line.

Efficiency varies with:

- Gamma energy
- Detector to sample geometry
- Matrix density
- Matrix composition

Experimental efficiency calibration in IUP laboratory

- In the IUP laboratory "standard geometries" using calibration sources with gamma-emitters mixture (wide range of energies)
 - liquids: Marinelli beakers, plastic bottles filled into defined heights ($1 \text{ g}\cdot\text{cm}^{-3}$)
 - solids: Marinelli beakers, plastic bottles filled into defined heights, filters and pressed pellets ($1, 1.5 \text{ or } 2 \text{ g}\cdot\text{cm}^{-3}$)

Experimental efficiency calibration in IUP laboratory

- In the IUP laboratory "standard geometries" using calibration sources with gamma-emitters mixture (wide range of energies)
 - liquids: Marinelli beakers, plastic bottles filled into defined heights ($1 \text{ g}\cdot\text{cm}^{-3}$)
 - solids: Marinelli beakers, plastic bottles filled into defined heights, filters and pressed pellets ($1, 1.5$ or $2 \text{ g}\cdot\text{cm}^{-3}$)

Sediment material:

- Different environments and different sources → enormous variability of compositions and densities.
- The matrix of the sample hardly comparable to the standard samples matrix → introducing additional error to efficiency calibration.

Experimental efficiency calibration in IUP laboratory

Density correction of samples with density different from "standard geometries":

- distributing the sample in a larger volume of inactive material (paraffin powder)
- preparing pressed pellets of geometries and densities close to those standard ones

Experimental efficiency calibration in IUP laboratory

Density correction of samples with density different from "standard geometries":

- distributing the sample in a larger volume of inactive material (paraffin powder)
- preparing pressed pellets of geometries and densities close to those standard ones

Suitable for:

- materials with higher activities (dilution does not matter)

Experimental efficiency calibration in IUP laboratory

Density correction of samples with density different from "standard geometries":

- distributing the sample in a larger volume of inactive material (paraffin powder)
- preparing pressed pellets of geometries and densities close to those standard ones


Suitable for:

- materials with higher activities (dilution does not matter)

Disadvantages for sediment samples:

- destructive method
- does not take into account changes of self absorption due to elemental composition of the sample matrix

Density and composition effects



Possible solutions for self absorption correction

- Transmission measurements for each sample (Cutshall, 1983; Hurtado et al., 2007) and applying corrections
- Mathematical efficiency calculation: full Monte Carlo calculations with MCNP, GEANT3, GEANT4, ...
- Specialized software, full MC: GESPECOR, DETEFF
- Specialized software, efficiency transfer: ETNA, EFFTRAN, ...
- Manufacturers' mathematical efficiency software: Ortec ANGLE, Canberra ISOCS, LabSOCS

Canberra detector characterization for LabSOCS

- LabSOCS = Laboratory SOurceless Calibration System
- Characterization of detectors:
 - model of the detector: physical dimensions, all internal and surrounding structures
 - modelling by MNCP - FEP efficiency response
 - model is compared to measurements using traceable sources in 5 different geometries - validation, adjustment of unknown parameters and dimensions (e.g., dead layers)
 - large number of efficiency datasets for point sources in vacuum at many locations around the detectors generated (20 energies between 10–7000 keV, 500 m diameter) and a calibration grid for a bare detector is created
 - attenuation due to internal structures of the detectors is calculated and attenuated efficiency grids at 20 energies → detector characterization file
- LabSOCS software calculates efficiency for voluminous samples by integrating the response over the volume of given source

LabSOCS Geometry Editor

A Bericht des Geometrie-Editors


Datum: Monday, June 02, 2008
Beschreibung: Det3_petri_soil_18mm_69.2_1.2
Kommentar:
Dateiname: c:\genie2k\cafiles-working-folder\gamma\an1\pitd
ana\testordner\geo\9501-4\geometry\det3_petri_soil_18mm_69.2_1.2.geo
Software: LabSOCS
Vorlage: SIMPLIFIED_BEAKER, Version: cone_or_cylinder
Detektor: S05007
Umgebung: Temperatur: 22 C, Druck: 760 mmHg, Rel. Luftfeuchte: 30%
Integrations: Konvergenz: 1.00%, MDRSN: 2*(4) CHSN: 2*(4)

# Geometrie-kompon.	Dimensionen (mm):						Material	D(g/cm ³)	R.Konz.
	d1	d2	d3	d4	d5	d6			
1 Becher	1.33	49.20	49.20	20.00			acrylic	1.17	
2 Probe	18.00						soil_1.0	1.20	
3 Adsorber1							none		
4 Adsorber2							none		
5 Source-Detektor									

Liste der Energien für die Erstellung der Effizienzkurve:

46.5 48.0 50.0 52.0 56.0 59.5 70.0 88.0 100.0 122.1
140.0 165.9 200.0 300.0 350.0 391.7 500.0 604.7 661.7 795.0
898.0 1173.2 1332.5 1460.8 1836.0

Simplified Beaker




Validation tests


A series of validation tests involving point and voluminous sources in different geometries.

- 1 a transmission test with point sources and different absorbers
- 2 point sources measured in different distances from the detector
- 3 a voluminous source measured in different distances from the detector
- 4 a soil sample measured in different geometries


Transmission test with ^{210}Pb and ^{137}Cs point sources and different absorbers



Point sources of ^{210}Pb and ^{137}Cs measured in different distances from the detector




A small comment to "small geometries"




Careful about detector's reference plain (detectors with carbon epoxy window tend to have concave top)

A voluminous source (^{137}Cs solution) measured in different distances from the detector



A soil sample with known ^{137}Cs and ^{40}K activity measured in different geometries



Summary and final remarks

- The claimed accuracy of LabSOCS efficiencies is from 4.3% for energies above 400 keV to 7.1% for 50–100 keV
- Mathematical calibration using LabSOCS: very useful tool, lots of freedom in sample characterization, BUT
- avoid "*blind confidence in factory-calibrated detector*" (Bossus et al., 2006; Gilmore, 2008)
- From Canberra recommended regular quality checks with test kit can only detect changes in detector parameters (growth of dead layer, etc.)
- Based on our experience, also testing LabSOCS efficiency calibration for voluminous samples for each detector is a good idea

Acknowledgments

Traveling grant: PIP (Physics International Postgraduate),
University of Bremen